## Leica Biosystems Differentiator products.

English.....	2
العربية (Arabic).....	6
简体中文 (Chinese Simplified).....	10
繁體中文 (Chinese Traditional).....	14
Dansk (Danish).....	18
Nederlands (Dutch).....	22
Français (French – Canada).....	26
Français (French – France).....	30
Deutsch (German).....	34
Italiano (Italian).....	38
日本語 (Japanese).....	41
한국어 (Korean).....	46
Norsk (Norwegian).....	50
Polski (Polish).....	54
Português (Portuguese – Brazil).....	58
Português (Portuguese – Portugal).....	62
Română (Romanian).....	66
Русский (Russian).....	70
Slovenščina (Slovenian).....	74
Español (Spanish – Central America).....	78
Español (Spanish – Spain).....	82
Svenska (Swedish).....	86
ภาษาไทย (Thai).....	90
Türkçe (Turkish).....	94
Tiếng Anh (Vietnamese).....	98

# Differentiators

**REF 3803590, 3803591, 3803595, 3803596, 3803598, 3803650, 3803650E, 3803651**

## Product Name

Leica Biosystems Differentiator products.

## Intended Use

### Detection/Measurement

The Leica Biosystems Differentiators do not detect or measure an analyte or marker. The Leica Biosystems Differentiators are intended to be used on frozen or paraffin embedded histology sections with a hematoxylin and eosin staining protocol. When used as recommended, the Differentiators remove excess of hematoxylin stain and define nuclei (in regressive staining method), remove non-specific hematoxylin staining from microscope glass slides caused by tissue adhesives, remove or minimize mucin staining providing transparency and definition of cytoplasmic stains.

### Product Function

The Leica Biosystems Differentiators are used in a H&E staining protocol to refine nuclei staining in regressive method and remove non-specific staining (e.g. glass slides or large, acidic proteins) in progressive and regressive staining methods. The H&E stained specimen, when interpreted by a trained professional, is utilized alongside other information such as the patient's medical history, physical condition, as well as results from other medical testing to render a medical diagnosis.

### Specific Information Provided

Leica Biosystems Differentiators are not intended for the detection, definition or differentiation of a specific disorder, condition or risk factor. The staining demonstrated with use of these products, when used as intended, provides trained professionals information which may define the physiological or pathological state of the tissue specimen.

### Automation

The Leica Biosystems Differentiators are not automated but can be used on automated staining platforms. Use on an automated platform should be validated at the point of use.

### Qualitative/Quantitative

The Leica Biosystems Differentiators are used with qualitative stains.

### Specimen Type

The Leica Biosystems Differentiators may be used with fixed or fresh histologic and cytologic specimens.

### Testing Population

The Leica Biosystems Differentiators are intended for use with any patient requiring evaluation of biopsy or resection tissue as for the assessment of a suspected pathology or disease.

### Intended User

The Leica Biosystems Differentiators are intended for use by qualified laboratory personnel and/or designee of the laboratory.

## In Vitro Diagnostic

The Leica Biosystems Differentiators are intended for *in vitro* diagnostics use only.

## Intended User

The Leica Biosystems Differentiators are intended for use by qualified laboratory personnel and/or designee.

## Test Principle

The Leica Biosystems Differentiators work by removing excess of hematoxylin stain from nuclei or any non-specific staining. The Differentiators are formulated with either strong or weak acid based on the staining methods used.

## Calibrators & Controls

The Leica Differentiators agents do not require the use of any calibrators or controls.

## Reagent Limitations

No reagent limitations are applicable to these products.

## Applicable Products

Product Code	Material Description
3803590	SelecTech Define Concentrate (500mL)
3803591	SelecTech Define Concentrate (4-500mL)
3803595	SelecTech Define MX-AQ Concentrate (500mL)
3803596	SelecTech Define MX-AQ Concentrate (4-500mL)
3803598	SelecTech Define MX-AQ RTU (Ready-To-Use) (1gal)
3803650	Surgipath Acid Alcohol 1% (1gal)
3803650E	Surgipath Acid Alcohol 0.5% (5L)
3803651	Surgipath Acid Alcohol 1% (4-1gal)

NOTE: Products listed here may not be available in all regions.

## Materials Not Included

The Leica Biosystems Differentiators are designed to be used as part of a Hematoxylin & Eosin (H&E) stain protocol which require the use of graded alcohols, xylene or xylene substitutes, hematoxylin, bluing agents, and eosin.

# Differentiators

**REF 3803590, 3803591, 3803595, 3803596, 3803598, 3803650, 3803650E, 3803651**

## Devices Required

Leica Biosystems Differentiators may be used on any open automated staining platform or with a manual staining method and should be validated at the point of use.

## Storage and Stability

The product shall be stable for 24 months postproduction when stored at ambient temperature.

Store reagents at room temperature (15-30°C) in a well-ventilated place.

**CAUTION:** Do not use after the expiration date.

## In Use Stability

User discretion should be utilized when determining in-use stability.

## Sterility

The Leica Biosystems Differentiators are not sterile products.

## Warnings/Precautions

This product and protocol(s) associated with the product, whether provided by Leica Biosystems in this instruction for use or developed by the user, shall be validated at the point of use by the user.

## Infectious Material Status

The Leica Biosystems Differentiators do not include any infectious material. However, specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions per facility guidelines.

## Special Facilities

The Leica Biosystems Differentiators should be used per facility guidelines.

## Specimen Handling

Specimens intended to be stained with H&E that includes Differentiators should be well fixed with neutral buffered formalin.

Following processing and paraffin embedding, section tissue at a standard thickness of (2 – 5µm).

## Preparation for Use

- **Define MX-aq RTU and Acid Alcohols (1% and 0.5%)** are ready-to-use formulas and therefore mixing is not required.
- **Define Concentrate and Define MX-aq Concentrate** must be diluted prior to use:
  - Define Concentrate: use one (1) part of concentrate, add nineteen (19) parts of 70% reagent alcohol to achieve working solution. To obtain a working solution, pour two measures of Concentrate into 500mL of 70% reagent alcohol and mix well. One measure equals fluid filled to first line in bottle cap.  
Note: Cap threads are not measure lines.
  - Define MX-aq Concentrate: use one (1) part of concentrate, add nineteen (19) parts of deionized or distilled water to achieve working solution. Alternatively, to obtain a working solution, pour two measures of concentrate into 500 mL of deionized or distilled water and mix well. One measure equals fluid filled to the first line in the bottle cap.  
Note: Cap threads are not measure lines.

Note - Inadequate washing after immersion in working solution of differentiator may affect desired staining intensity.

## Protocol Set-up:

1. Use in accordance with established hematoxylin and eosin staining procedures.
2. Following hematoxylin staining, rinse slides with running tap water until all free hematoxylin is removed.
3. Immerse slides in a working differentiator solution for a time period as necessary to remove background hematoxylin staining from the slide or tissue.  
Immersion times in a differentiator depend on the acid type:
  - Immerse in Define or Define MX-aq, both containing weak organic acid, for thirty (30) seconds to ninety (90) seconds. Example of progressive staining protocol with regressive element in shown in Image 1.
  - Immerse in 1% or 0.5% Acid Alcohols for three (3) to ten (10) seconds. This is differentiating when using regressive staining protocol (Image 2). Example of regressive staining protocol in shown in Image 2.
4. Rinse in running tap water for thirty (30) seconds to one (1) minute.
5. Continue with hematoxylin and eosin staining protocol.

# Differentiators

REF 3803590, 3803591, 3803595, 3803596, 3803598, 3803650, 3803650E, 3803651

**Table 1. Example of progressive H&E staining protocol with regressive element, utilizing weak organic acid (Define).**

Steps	Action	Chemical	Time (mm:ss)
1	Deparaffinize	Xylene	3:00
2	Deparaffinize	Xylene	3:00
3	Deparaffinize	Xylene	3:00
4	Hydration	100% Alcohol	2:00
5	Hydration	100% Alcohol	1:00
6	Hydration	100% Alcohol	1:00
7	Hydration	80% or 95% Alcohol	1:00
8	Hydration	Water Wash	1:00
9	Stain	Progressive Hematoxylin	1:00 to 5:00
10	Wash	Water Wash	3:00
11	Differentiation	Differentiator	0:30 to 1:30
12	Wash	Water Wash	1:00
13	Bluing	Bluing Buffer	0:30 to 1:00
14	Wash	Water Wash	2:00
15	Dehydration	80% to 95% Alcohol	1:00
16	Counterstaining	Eosin	0:30 to 1:30
17	Wash	Water Wash	2:00
18	Dehydration	95% to 100% Alcohol	1:00
19	Dehydration	100% Alcohol	1:00
20	Dehydration	100% Alcohol	1:00
21	Clearing	Xylene	2:00
22	Clearing	Xylene	2:00
23	Clearing	Xylene	2:00

**Table 2. Example of regressive H&E staining protocol utilizing strong inorganic acid (Acid Alcohol).**

Steps	Action	Chemical	Time (mm:ss)
1	Deparaffinize	Xylene	3:00
2	Deparaffinize	Xylene	3:00
3	Deparaffinize	Xylene	3:00
4	Hydration	100% Alcohol	2:00
5	Hydration	100% Alcohol	1:00
6	Hydration	100% Alcohol	1:00
7	Hydration	80% or 95% Alcohol	1:00
8	Hydration	Water Wash	1:00
9	Stain	Progressive Hematoxylin	1:00 to 5:00
10	Wash	Water Wash	3:00
11	Differentiation	Differentiator	0:03 to 0:10
12	Wash	Water Wash	1:00
13	Bluing	Bluing Buffer	0:30 to 1:00
14	Wash	Water Wash	2:00
15	Dehydration	80% to 95% Alcohol	1:00
16	Counterstaining	Eosin	0:30 to 1:30
17	Wash	Water Wash	2:00
18	Dehydration	95% to 100% Alcohol	1:00

# Differentiators

**REF 3803590, 3803591, 3803595, 3803596, 3803598, 3803650, 3803650E, 3803651**

19	Dehydration	100% Alcohol	1:00
20	Dehydration	100% Alcohol	1:00
21	Clearing	Xylene	2:00
22	Clearing	Xylene	2:00
23	Clearing	Xylene	2:00

### Readiness for Use

Once concentrated version of a differentiator is properly diluted, or if ready-to-use formula is used, pour all the reagent into the reagent vessel. Place the reagent vessel back into the respective station.

### Quality Control

A routine quality control slide(s) containing tissue fixed and processed in a similar manner to the test specimens should be performed prior to routine use to ensure reagents are performing as intended.

### Expected Results

By following the instruction for use, differentiators shall remove excessive or background/non-specific staining from hematoxylin stained specimens such that nuclear detail shall be visualized microscopically.

### Analytical Performance

The Leica Biosystems Differentiators are not used to detect a specific analyte or marker. These products are used in conjunction with other products in a Hematoxylin & Eosin staining protocol system to stain cell nuclei blue and connective tissue, cytoplasm, muscle and erythrocytes various shades of orange, pink and red. Analytical parameters such as analytical sensitivity, analytical specificity, trueness (bias), precision (repeatability and reproducibility), accuracy (resulting from trueness and precision), limits of detection and quantitation, measuring range, linearity, cut-off, including determination of appropriate criteria for specimen collection and handling and control of known relevant endogenous and exogenous interference, cross-reactions do not apply to the performance of this system.

### Clinical Performance

The Leica Biosystems Differentiators are not intended for use as a means of detecting a specific disease or pathological process or state. Clinical performance indices such as diagnostic sensitivity, diagnostic specificity, positive predictive value, negative predictive value, likelihood ratio as well as expected values in normal and affected populations do not apply to the use of the Leica Biosystems Differentiators in a clinical setting.

### Disposal

Differentiators should be disposed in accordance with local governing regulations.

Leica Biosystems Richmond, Inc.  
5205 Route 12  
Richmond, IL 60071  
USA  
(1-844-534-2262)  
LeicaBiosystems.com

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