BIO SYSTEMS

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BOND™ Aspirating Probe Cleaning System

Catalog No: CS9100

Intended Use

The BOND Aspirating Probe Cleaning System is intended for automated cleaning of aspirating probes on BOND instruments, using the aspirating probe cleaning protocol installed in BOND systems (includes Leica BOND-MAX system and Leica BOND-III system).

Summary and Explanation

Routine use of the BOND Aspirating Probe Cleaning System keeps BOND system aspirating probes free of DAB residue.

The BOND Aspirating Probe Cleaning System is a tray of three cleaning reagents that is loaded onto the BOND system like any other set of reagents. A 20-minute predefined cleaning protocol mixes reagents from the system in the BOND mixing station then uses a sequence of aspirations and dispenses in the station to thoroughly clean the bottom 2–3 cm of the probe. As a final step the probe is rinsed with BOND Wash to remove any cleaning reagent residue.

Use of the cleaning system does not remove the need for regular, manual cleaning of the probe. Leica Biosystems recommends that users wipe the external surface of the probe with a soft cloth dipped in 70% alcohol once a week.

Reagents Provided

Reagents sufficient for 15 cleans:

1. Cleaning Component A (3.75 mL) proprietary cleaning reagent
2. Cleaning Component B (3.75 mL) <5% sulphuric acid
3. Cleaning Component C (3.75 mL) proprietary cleaning reagent

Dilution and Mixing

The BOND Aspirating Probe Cleaning System is optimized for use on the BOND system. No manual reconstitution, mixing, dilution, or titration is required. BOND mixes components to the correct concentrations on the instrument when cleaning is initiated.

Materials Required but Not Provided

No additional materials are required to use the BOND Aspirating Probe Cleaning System.

Storage and Stability

Store at 2–8 °C. The product is stable under these conditions up to the expiry date indicated on the container and tray handle labels.

There are no obvious signs to indicate instability of this product.

Return to 2–8 °C immediately after use.

Storage conditions other than those specified above must be verified by the user.

Precautions

- This product is intended for *in vitro* diagnostic use.

CLEANING COMPONENT A

Contains Potassium Permanganate (<1%).
GHS09: Environmental.

H411: Toxic to aquatic life with long lasting effects.

P273: Avoid release to the environment.
P391: Collect spillage.

CLEANING COMPONENT B

Contains Sulphuric Acid (<10%).
GHS07: Exclamation mark.
Signal words: Warning.

H315: Causes skin irritation.
H319: Causes serious eye irritation.

P264: Wash hands thoroughly after handling.
P280: Wear protective gloves/protective clothing/eye protection/face protection.
P302+352: IF ON SKIN: Wash with plenty of soap and water.
P332+313: If skin irritation occurs: Get medical advice/attention.
P362: Take off contaminated clothing and wash before reuse.
P305+351+338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P337+313: If eye irritation persists: Get medical advice/attention.

CLEANING COMPONENT C

Contains Oxalic Acid (<10%).
GHS07: Exclamation mark.
Signal words: Warning.

H302 + H312: Harmful if swallowed or in contact with skin.

P264: Wash hands thoroughly after handling.
P270: Do not eat, drink or smoke when using this product.
P280: Wear protective gloves/protective clothing/eye protection/face protection.
P301 + P312: IF SWALLOWED: Call a POISON CENTRE or doctor if you feel unwell.
P330: Rinse mouth.
P302 + P352: IF ON SKIN: Wash with plenty of soap and water.
P312: Call a POISON CENTRE or doctor if you feel unwell.
P363: Wash contaminated clothing before reuse.
P501: Dispose of contents/container to hazardous or special waste collection point.

- If a spill occurs on-board the BOND system, it should be cleaned up as soon as possible using an appropriate absorbent material.
- To obtain a copy of the Material Safety Data Sheet (MSDS) contact your local distributor or Leica Biosystems' regional office, or visit the Leica Biosystems, Biosystems Division Web site: www.LeicaBiosystems.com
- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions¹. Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.
- Consult federal, state and local regulations for disposal of any potentially toxic components.

Instructions for Use

Follow the instructions below to clean the aspirating probe with the BOND Aspirating Probe Cleaning System.

The cleaning protocol takes approximately 20 minutes to run.

1. Ensure that the Processing Module is idle, with no trays loaded, scheduled, or running.
2. Insert the BOND Aspirating Probe Cleaning System into the reagent rack on the Processing Module.
3. Initiate the clean by clicking "Clean aspirating probe" in the submenu for the particular Processing Module in the Maintenance menu.
4. Click Yes to start the clean when prompted.
5. The cleaning protocol begins, indicated by the cleaning icon in the Processing Module tab.
6. Wait until notified that the clean is finished.
7. Remove the BOND Aspirating Probe Cleaning System from the reagent rack.
8. Click OK in the cleaning completion notification dialog to continue normal operation.

References

1. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.

Date of Issue

10 September 2018

BOND™ Aspirating Probe Cleaning System

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Leica Biosystems Newcastle Ltd
Leica Biosystems Inc
Pty Ltd