## BOND Polymer Refine Detection

### Catalog No:DS9800

Leica Biosystems Newcastle Ltd  
Balliol Business Park West  
Benton Lane  
Newcastle Upon Tyne NE12 8EW  
United Kingdom  
☎ +44 191 215 4242

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## Catalog No: DS9800

### Intended Use

This detection system is for *in vitro* diagnostic use.

BOND Polymer Refine Detection is a biotin-free, polymeric horseradish peroxidase (HRP)-linker antibody conjugate system for the detection of tissue-bound mouse and rabbit IgG and some mouse IgM primary antibodies. It is intended for staining sections of formalin-fixed, paraffin-embedded tissue on the BOND automated system (includes Leica BOND-MAX system and Leica BOND-III system).

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls. They should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

The BOND Polymer Refine Detection Kit must be used with laboratory best practice in the use of tissue controls. For assurance, laboratories should stain each patient sample in conjunction with positive, negative, and other tissue specific controls as needed.

### Summary and Explanation

Immunohistochemical techniques can be used to demonstrate the presence of antigens in tissue and cells (see "Using BOND Reagents" in your BOND user documentation).

BOND Polymer Refine Detection utilizes a novel controlled polymerization technology to prepare polymeric HRP-linker antibody conjugates. The detection system avoids the use of streptavidin and biotin, and therefore eliminates non-specific staining as a result of endogenous biotin.

BOND Polymer Refine Detection works as follows:

- The specimen is incubated with hydrogen peroxide to quench endogenous peroxidase activity.
- A user-supplied specific primary antibody is applied.
- Post Primary IgG linker reagent localizes mouse antibodies.
- Poly-HRP IgG reagent localizes rabbit antibodies.
- The substrate chromogen, 3,3'-Diaminobenzidine tetrahydrochloride hydrate (DAB), visualizes the complex via a brown precipitate.
- Hematoxylin (blue) counterstaining allows the visualization of cell nuclei.

Using BOND Polymer Refine Detection in combination with the BOND automated system reduces the possibility of human error and inherent variability resulting from individual reagent dilution, manual pipetting and reagent application.

### Reagents Provided

Reagents sufficient for 200–300 tests

1. Peroxide Block (30 mL) 3–4% (v/v) Hydrogen peroxide.
2. Post Primary (30 mL) Rabbit anti mouse IgG (<10 µg/mL) in 10% (v/v) animal serum in tris-buffered saline/0.1% ProClin™ 950.
3. Polymer (30 mL) Anti-rabbit Poly-HRP-IgG (<25 µg/mL) containing 10% (v/v) animal serum in tris-buffered saline/0.1% ProClin™ 950.
4. DAB Part 1 (2.4 mL) 66 mM 3,3'-Diaminobenzidine tetrahydrochloride hydrate, in a stabilizer solution.
5. DAB Part B (30 mL) ≤0.1% (v/v) Hydrogen Peroxide in a stabilizer solution.
6. DAB Part B (30 mL) ≤0.1% (v/v) Hydrogen Peroxide in a stabilizer solution.
7. Hematoxylin (30 mL) <0.1% Hematoxylin.

### Dilution and Mixing

BOND Polymer Refine Detection is optimized for use on the BOND system. Reconstitution, mixing, dilution, or titration of these reagents is not required.

### Materials Required But Not Provided

Refer to "Using BOND Reagents" in your BOND user documentation for a complete list of materials required for specimen treatment and immunohistochemical staining using the BOND system.

### Storage and Stability

Store at 2–8 °C. Do not freeze. Do not use after the expiration date indicated on the tray handle label. Return to 2–8 °C immediately after use.

There are no obvious signs to indicate instability of this product, therefore positive and negative controls should be run simultaneously with unknown specimens (refer to "Quality Control" in the "Using BOND Reagents" section of your BOND user documentation).

If unexpected staining is observed that cannot be explained by variations in laboratory procedures, and a problem with the detection system is suspected, contact your local distributor or the regional office of Leica Biosystems immediately.

Storage conditions other than those specified above must be verified by the user<sup>1</sup>.

## Precautions

- Restricted to professional users.
- This detection system is intended for *in vitro* diagnostic use.
- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems' web site, [www.LeicaBiosystems.com](http://www.LeicaBiosystems.com)
- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions<sup>3</sup>. Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.
- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user<sup>1</sup>.
- Do not mix reagents from different detection systems.

## Instructions for Use

BOND Polymer Refine Detection was developed for use on the BOND automated system. Operating parameters for application of the detection system reagents on the BOND Processing Module have been optimized at Leica Biosystems. These can be displayed by following the instructions in your BOND user documentation.

## Product Specific Limitations

BOND Polymer Refine Detection has been optimized at Leica Biosystems for use with BOND ancillary reagents. Laboratories may use their own primary antibodies provided they have been diluted to an appropriate concentration with BOND Primary Antibody Diluent (Catalog No. AR9352). Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances.

The appropriate concentration of user's own primary antibodies may vary, due to variation in tissue fixation and the effectiveness of antigen enhancement, and must be determined empirically. Negative reagent controls should be used when optimizing retrieval conditions and primary antibody concentrations.

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls. They should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

## Troubleshooting

Refer to reference 4 for remedial action.

If the patient result does not correspond to the expected results with the use of controls, the test should be repeated.

If the staining result is not as expected, and you wish to troubleshoot performance of the instrument and detection system independently, your local Leica representative can provide specific protocols. The detection kit must be used in conformance to the package instructions and within the shelf life indicated on the product itself.

## Further Information

Further information on immunostaining with BOND reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels, and General Limitations can be found in "Using BOND Reagents" in your BOND user documentation.

## Bibliography

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. List of substances which may be candidates for further scientific review and possible identification, classification and regulation as potential occupational carcinogens. Fed Reg 1980; 45:157.
3. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order Code M29-P.
4. JD Bancroft and A Stevens. Theory and Practice of Histological Techniques. 4<sup>th</sup> Edition. Churchill Livingstone, New York. 1996.

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## Bibliographie

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. List of substances which may be candidates for further scientific review and possible identification, classification and regulation as potential occupational carcinogens. Fed Reg 1980; 45:157.
3. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
4. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.

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## Bibliografia

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. List of substances which may be candidates for further scientific review and possible identification, classification and regulation as potential occupational carcinogens. Fed Reg 1980; 45:157.
3. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order Code M29-P.
4. JD Bancroft and A Stevens. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.

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1. Clinical laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. List of substances which may be candidates for further scientific review and possible identification, classification and regulation as potential occupational carcinogens. Fed Reg 1980; 45:157.
3. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order Code M29-P.
4. JD Bancroft and A Stevens. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.

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## Bibliografi

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. List of substances which may be candidates for further scientific review and possible identification, classification and regulation as potential occupational carcinogens. Fed Reg 1980; 45:157.
3. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order Code M29-P.
4. JD Bancroft and A Stevens. Theory and Practice of Histological Techniques. 4<sup>th</sup> Edition Churchill Livingstone, New York. 1996.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. List of substances which may be candidates for further scientific review and possible identification, classification and regulation as potential occupational carcinogens. Fed Reg 1980; 45:157.
3. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
4. JD Bancroft and A Stevens. Theory and Practice of Histological Techniques. 4<sup>th</sup> Edition. Churchill Livingstone, New York. 1996.

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Leica Biosystems Inc
Leica Biosystems Melbourne Pty Ltd