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# BOND™ Ready-To-Use Primary Antibody CD56 (CD564)

Catalog No: PA0191

## Intended Use

This reagent is for *in vitro* diagnostic use.

CD56 (CD564) monoclonal antibody is intended to be used for the qualitative identification by light microscopy of human CD56 molecules (neural cell adhesion molecules) in formalin-fixed, paraffin-embedded tissue by immunohistochemical staining using the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

## Summary and Explanation

Immunohistochemical techniques can be used to demonstrate the presence of antigens in tissue and cells (see "Using BOND Reagents" in your BOND user documentation). CD56 (CD564) primary antibody is a ready to use product that has been specifically optimized for use with BOND Polymer Refine Detection. The demonstration of human CD56 molecules (neural cell adhesion molecules) is achieved by first allowing the binding of CD56 (CD564) to the section, and then visualizing this binding using the reagents provided in the detection system. The use of these products, in combination with the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system), reduces the possibility of human error and inherent variability resulting from individual reagent dilution, manual pipetting and reagent application.

## Reagents Provided

CD56 (CD564) is a mouse anti-human monoclonal antibody produced as a tissue culture supernatant, and supplied in Tris buffered saline with carrier protein, containing 0.35 % ProClin™ 950 as a preservative.

Total volume = 7 mL.

## Clone

CD564

## Immunogen

Prokaryotic recombinant fusion protein corresponding to the external domain of the human CD56 molecule.

## Specificity

Human CD56 molecule (neural cell adhesion molecule).

## Ig Class

IgG2b

## Total Protein Concentration

Approx 10 mg/mL.

## Antibody Concentration

Greater than or equal to 0.24 mg/L as determined by ELISA.

## Dilution and Mixing

CD56 (CD564) primary antibody is optimally diluted for use on the BOND system (includes Leica BOND-MAX system and Leica BOND-III system). Reconstitution, mixing, dilution or titration of this reagent is not required.

## Materials Required But Not Provided

Refer to "Using BOND Reagents" in your BOND user documentation for a complete list of materials required for specimen treatment and immunohistochemical staining using the BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

## Storage and Stability

Store at 2–8 °C. Do not use after the expiration date indicated on the container label.

The signs indicating contamination and/or instability of CD56 (CD564) are: turbidity of the solution, odor development, and presence of precipitate.

Return to 2–8 °C immediately after use.

Storage conditions other than those specified above must be verified by the user<sup>1</sup>.

## Precautions

- This product is intended for *in vitro* diagnostic use.
- The concentration of ProClin™ 950 is 0.35 %. It contains the active ingredient 2-methyl-4-isothiazolin-3-one, and may cause irritation to the skin, eyes, mucous membranes and upper respiratory tract. Wear disposable gloves when handling reagents.
- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems' Web site, [www.LeicaBiosystems.com](http://www.LeicaBiosystems.com)

- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions<sup>2</sup>. Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.
- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Retrieval, incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user.

## Instructions for Use

CD56 (CD564) primary antibody was developed for use on the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system) in combination with BOND Polymer Refine Detection. The recommended staining protocol for CD56 (CD564) primary antibody is IHC Protocol F. Heat induced epitope retrieval is recommended using BOND Epitope Retrieval Solution 1 for 20 minutes.

## Results Expected

### Normal Tissues

Clone CD564 detected the CD56 (neural cell adhesion molecule) antigen on the membrane of NK cells, a subset of activated T cells and neuroectodermal cells. (Total number of normal cases evaluated = 79).

### Tumor Tissues

Clone CD564 stained 6/6 NK/T cell lymphomas, 3/115 diffuse large B-cell lymphomas, 3/9 plasma cell myelomas, 2/8 T-cell lymphomas, 2/2 neuroblastomas, 1/3 Burkitt's lymphomas, 1/1 primitive B/T cell acute lymphoblastic lymphoma, 3/3 carcinoid tumors, 3/3 medullary carcinomas of the thyroid, 2/6 lung tumors (including 1/2 small cell carcinomas, 1/1 non-small cell carcinoma, 0/1 adenocarcinoma, 0/1 squamous cell carcinoma and 0/1 large cell carcinoma), 2/2 Wilms' tumors, 1/1 pheochromocytoma, 1/1 anaplastic astrocytoma of the brain, 1/1 choroid plexus papilloma of the brain, 1/1 atypical carcinoid of the thymus, 1/1 soft tissue ganglioneuroma, 1/2 kidney renal cell carcinomas, 1/2 seminomas, 1/4 liver carcinomas (including 1/1 liver metastatic carcinoma, 0/1 cholangiocarcinoma and 0/2 hepatocellular carcinomas) and 1/4 ovarian carcinomas (including 1/1 serous cystadenocarcinoma, 0/1 malignant germ cell tumor, 0/1 clear cell carcinoma and 0/1 mucinous cystadenocarcinoma). No staining was observed in Hodgkin's disease (0/20), chronic lymphocytic lymphomas (0/12), follicular lymphomas (0/11), T-cell anaplastic large cell lymphomas (0/7), mantle cell lymphomas (0/7), angioimmunoblastic T-cell lymphomas (0/4), acute myeloid leukemias (0/3), follicular B-cell lymphomas (0/2), a B-cell acute lymphoblastic lymphoma (0/1), a peripheral T-cell lymphoma (0/1), a marginal zone lymphoma (0/1), a myeloid leukemia (0/1), thyroid tumors (0/4), esophageal tumors (0/2), breast tumors (0/2), metastatic tumors of unknown origin (0/2), tumors of the tongue (0/2), tumors of the stomach (0/2), tumors of the colon (0/2), rectal tumors (0/2), cervical tumors (0/2), tumors of the skin (0/2), a soft tissue fibromatosis (0/1), a tumor of the larynx (0/1), a plasmacytoma (0/1), a neuroectodermal tumor (0/1) or a Ewing's sarcoma (0/1). (Total number of tumor cases evaluated = 272).

**CD56 (CD564) is recommended for the detection of human CD56 protein in normal and neoplastic tissues, as an adjunct to conventional histopathology using non-immunologic histochemical stains.**

## Product Specific Limitations

CD56 (CD564) has been optimized at Leica Biosystems for use with BOND Polymer Refine Detection and BOND ancillary reagents. Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances. The protocol times may vary, due to variation in tissue fixation and the effectiveness of antigen enhancement, and must be determined empirically. Negative reagent controls should be used when optimizing retrieval conditions and protocol times.

## Troubleshooting

Refer to reference 3 for remedial action.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

## Further Information

Further information on immunostaining with BOND reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels, and General Limitations can be found in "Using BOND Reagents" in your BOND user documentation.

## Bibliography

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Nakamoto Y, Mizukoshi E, Tsuji H, et al. Combined therapy of transcatheter hepatic arterial embolization with intratumoral dendritic cell infusion for hepatocellular carcinoma: clinical safety. Clinical and Experimental Immunology. 2007; 147:296-305.
5. Wicherek L, Dutsch-Wicherek M, Galazka K, et al. Comparison of RCAS1 and metallothionein expression and the presence and activity of immune cells in human ovarian and abdominal wall endometriomas. Reproductive Biology and Endocrinology. 2006; 4:41.
6. Wicherek I, Galazka K, Popiela TJ, et al. Metallothionein expression and infiltration of cytotoxic lymphocytes in uterine and tubal implantation sites. Journal of Reproductive Immunology. 2006; 70:119-131.

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## Bibliografie

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Nakamoto Y, Mizukoshi E, Tsuji H, et al. Combined therapy of transcatheter hepatic arterial embolization with intratumoral dendritic cell infusion for hepatocellular carcinoma: clinical safety. Clinical and Experimental Immunology. 2007; 147:296-305.
5. Wicherek L, Dutsch-Wicherek M, Galazka K, et al. Comparison of RCAS1 and metallothionein expression and the presence and activity of immune cells in human ovarian and abdominal wall endometriomas. Reproductive Biology and Endocrinology. 2006; 4:41.
6. Wicherek I, Galazka K, Popiela T.J, et al. Metallothionein expression and infiltration of cytotoxic lymphocytes in uterine and tubal implantation sites. Journal of Reproductive Immunology. 2006; 70:119-131.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Nakamoto Y, Mizukoshi E, Tsuji H, et al. Combined therapy of transcatheter hepatic arterial embolization with intratumoral dendritic cell infusion for hepatocellular carcinoma: clinical safety. Clinical and Experimental Immunology. 2007; 147:296-305.
5. Wicherek L, Dutsch-Wicherek M, Galazka K, et al. Comparison of RCAS1 and metallothionein expression and the presence and activity of immune cells in human ovarian and abdominal wall endometriomas. Reproductive Biology and Endocrinology. 2006; 4:41.
6. Wicherek I, Galazka K, Popielka T, et al. Metallothionein expression and infiltration of cytotoxic lymphocytes in uterine and tubal implantation sites. Journal of Reproductive Immunology. 2006; 70:119-131.

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## Bibliografi

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Nakamoto Y, Mizukoshi E, Tsuji H, et al. Combined therapy of transcatheter hepatic arterial embolization with intratumoral dendritic cell infusion for hepatocellular carcinoma: clinical safety. Clinical and Experimental Immunology. 2007; 147:296-305.
5. Wicherek L, Dutsch-Wicherek M, Galazka K, et al. Comparison of RCAS1 and metallothionein expression and the presence and activity of immune cells in human ovarian and abdominal wall endometriomas. Reproductive Biology and Endocrinology. 2006; 4:4-11.
6. Wicherek I, Galazka K, Popiela TJ, et al. Metallothionein expression and infiltration of cytotoxic lymphocytes in uterine and tubal implantation sites. Journal of Reproductive Immunology. 2006; 70:119-131.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Nakamoto Y, Mizukoshi E, Tsuji H, et al. Combined therapy of transcatheter hepatic arterial embolization with intratumoral dendritic cell infusion for hepatocellular carcinoma: clinical safety. Clinical and Experimental Immunology. 2007; 147:296-305.
5. Wicherek L, Dutsch-Wicherek M, Galazka K, et al. Comparison of RCAS1 and metallothionein expression and the presence and activity of immune cells in human ovarian and abdominal wall endometriomas. Reproductive Biology and Endocrinology. 2006; 4:4-1.
6. Wicherek I, Galazka K, Popiela T, et al. Metallothionein expression and infiltration of cytotoxic lymphocytes in uterine and tubal implantation sites. Journal of Reproductive Immunology. 2006; 70:119-131.

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## Kaynakça

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Nakamoto Y, Mizukoshi E, Tsuji H, et al. Combined therapy of transcatheter hepatic arterial embolization with intratumoral dendritic cell infusion for hepatocellular carcinoma: clinical safety. Clinical and Experimental Immunology. 2007; 147:296-305.
5. Wicherek L, Dutsch-Wicherek M, Galazka K, et al. Comparison of RCA1 and metallothionein expression and the presence and activity of immune cells in human ovarian and abdominal wall endometriomas. Reproductive Biology and Endocrinology. 2006; 4:41.
6. Wicherek I, Galazka K, Popiela TJ, et al. Metallothionein expression and infiltration of cytotoxic lymphocytes in uterine and tubal implantation sites. Journal of Reproductive Immunology. 2006; 70:119-131.

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## Bibliografie

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Nakamoto Y, Mizukoshi E, Tsuji H, et al. Combined therapy of transcatheter hepatic arterial embolization with intratumoral dendritic cell infusion for hepatocellular carcinoma: clinical safety. Clinical and Experimental Immunology. 2007; 147:296-305.
5. Wicherek L, Dutsch-Wicherek M, Galazka K, et al. Comparison of RCAS1 and metallothionein expression and the presence and activity of immune cells in human ovarian and abdominal wall endometriosis. Reproductive Biology and Endocrinology. 2006; 4:41.
6. Wicherek I, Galazka K, Popiela TJ, et al. Metallothionein expression and infiltration of cytotoxic lymphocytes in uterine and tubal implantation sites. Journal of Reproductive Immunology. 2006; 70:119-131.

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## Literatúra

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Nakamoto Y, Mizukoshi E, Tsuji H, et al. Combined therapy of transcatheter hepatic arterial embolization with intratumoral dendritic cell infusion for hepatocellular carcinoma: clinical safety. Clinical and Experimental Immunology. 2007; 147:296-305.
5. Wicherek L, Dutsch-Wicherek M, Galazka K, et al. Comparison of RCAS1 and metallothionein expression and the presence and activity of immune cells in human ovarian and abdominal wall endometriosis. Reproductive Biology and Endocrinology. 2006; 4:4.1.
6. Wicherek I, Galazka K, Popiela T.J, et al. Metallothionein expression and infiltration of cytotoxic lymphocytes in uterine and tubal implantation sites. Journal of Reproductive Immunology. 2006; 70:119-131.

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Leica Biosystems Newcastle Ltd
Leica Biosystems Inc
Pty Ltd