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BOND™ Ready-To-Use Primary Antibody CD15 (MMA)

Catalog No: PA0473

Intended Use

This reagent is for *in vitro* diagnostic use.

CD15 (MMA) monoclonal antibody is intended to be used for the qualitative identification by light microscopy of human CD15 protein in formalin-fixed, paraffin-embedded tissue by immunohistochemical staining using the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Summary and Explanation

Immunohistochemical techniques can be used to demonstrate the presence of antigens in tissue and cells (see "Using BOND Reagents" in your BOND user documentation). CD15 (MMA) primary antibody is a ready to use product that has been specifically optimized for use with BOND Polymer Refine Detection. The demonstration of human CD15 protein is achieved by first allowing the binding of CD15 (MMA) to the section, and then visualizing this binding using the reagents provided in the detection system. The use of these products, in combination with the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system), reduces the possibility of human error and inherent variability resulting from individual reagent dilution, manual pipetting and reagent application.

Reagents Provided

CD15 (MMA) is a mouse anti-human monoclonal antibody produced as a tissue culture supernatant, and supplied in Tris buffered saline with carrier protein, containing 0.35 % ProClin™ 950 as a preservative.

Total volume = 7 mL.

Clone

MMA.

Immunogen

BALB/C mice injected with U937 histiocytic cell line.

Specificity

Human CD15 antigen.

Ig Class

IgM.

Total Protein Concentration

Approx 10 mg/mL.

Antibody Concentration

Not applicable.

Dilution and Mixing

CD15 (MMA) primary antibody is optimally diluted for use on the BOND system (includes Leica BOND-MAX system and Leica BOND-III system). Reconstitution, mixing, dilution or titration of this reagent is not required.

Materials Required But Not Provided

Refer to "Using BOND Reagents" in your BOND user documentation for a complete list of materials required for specimen treatment and immunohistochemical staining using the BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

Storage and Stability

Store at 2–8 °C. Do not use after the expiration date indicated on the container label.

The signs indicating contamination and/or instability of CD15 (MMA) are: turbidity of the solution, odor development, and presence of precipitate.

Return to 2–8 °C immediately after use.

Storage conditions other than those specified above must be verified by the user¹.

Precautions

- This product is intended for *in vitro* diagnostic use.
- The concentration of ProClin™ 950 is 0.35 %. It contains the active ingredient 2-methyl-4-isothiazolin-3-one, and may cause irritation to the skin, eyes, mucous membranes and upper respiratory tract. Wear disposable gloves when handling reagents.
- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems' Web site, www.LeicaBiosystems.com

- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions². Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.
- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Retrieval, incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user.

Instructions for Use

CD15 (MMA) primary antibody was developed for use on the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system) in combination with BOND Polymer Refine Detection. The recommended staining protocol for CD15 (MMA) primary antibody is IHC Protocol F. Heat induced epitope retrieval is recommended using BOND Epitope Retrieval Solution 1 for 20 minutes.

Results Expected

Normal Tissues

Clone MMA detected the CD15 protein on the cell membrane and in the cytoplasm of epithelial cells and granulocytes in a variety of tissues evaluated. Staining was also observed in neuropil of the cerebrum and cerebellum, Leydig cells in testis and Hassall's corpuscles in thymus (Total number of normal cases evaluated = 123).

Tumor Tissues

Clone MMA stained 36/195 hematological malignancies (including 34/40 Hodgkin's lymphomas, 1/107 diffuse large B-cell lymphomas, 1/1 peripheral T-cell lymphoma, 0/11 chronic lymphocytic lymphomas, 0/11 follicular lymphomas, 0/7 T-cell anaplastic large cell lymphomas, 0/6 mantle cell lymphomas, 0/4 angioimmunoblastic T-cell lymphomas, 0/3 T/NK cell lymphomas, 0/1 B-cell acute lymphoblastic lymphoma, 0/1 primitive B/T cell acute lymphoblastic lymphoma, 0/1 T-cell lymphoma, 0/1 marginal zone lymphoma and 0/1 non-Hodgkin's B-cell lymphoma). Clone MMA also stained 8/9 bowel tumors, 3/4 lung tumors, 2/3 adenocarcinomas of the stomach, 2/2 clear cell carcinomas of the kidney, 2/2 ovarian tumors, 2/2 endometrial tumors, 1/5 brain tumors, 1/5 metastatic tumors, 1/3 squamous cell carcinomas of the esophagus, 1/2 tumors of the adrenal gland and 1/2 tumors of the head and neck. Except for infiltrating granulocytes, no staining was detected in a variety of additional abnormal tissues evaluated, including breast tumors (0/5), liver tumors (0/5), tumors of the thyroid (0/5), bladder tumors (0/2), prostatic tumors (0/2), melanomas (0/2), seminomas (0/2), cervical tumors (0/2), a tongue tumor (0/1), a pancreatic tumor (0/1), a prostatic hyperplasia (0/1), a tumor of the salivary gland (0/1), a skin tumor (0/1) and a bone tumor (0/1). (Total number of abnormal cases evaluated = 265).

CD15 (MMA) is recommended for the detection of CD15 protein in normal and neoplastic tissues, as an adjunct to conventional histopathology using non-immunologic histochemical stains.

Product Specific Limitations

CD15 (MMA) has been optimized at Leica Biosystems for use with BOND Polymer Refine Detection and BOND ancillary reagents. Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances. The protocol times may vary, due to variation in tissue fixation and the effectiveness of antigen enhancement, and must be determined empirically. Negative reagent controls should be used when optimizing retrieval conditions and protocol times.

Troubleshooting

Refer to reference 3 for remedial action.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

Further Information

Further information on immunostaining with BOND reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels, and General Limitations can be found in "Using BOND Reagents" in your BOND user documentation.

Bibliography

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Mao X, Zhang X, Xue X, et al. Brain tumor stem-like cells identified by neural stem cell marker CD15. *Translational Oncology*. 2009; 2(4): 247-257.
5. Wong A, Lopategui J, Clancy S, et al. Anaplastic large cell lymphoma associated with a breast implant capsule: a case report and review of the literature. *The American Journal of Surgical Pathology*. 2008; 32(8):1265-1268.
6. Pellegrini W, Bresciani G, De Zorzi A, et al. MMA monoclonal antibody is a superior anti-CD15 reagent for the diagnosis of classical Hodgkin's lymphoma? *Haematologica*. 2007; 92(5):708-709.
7. Vassallo J, Lamant L, Brugieres L, et al. ALK positive anaplastic large cell lymphoma mimicking nodular sclerosis Hodgkin's lymphoma: Report of 10 cases. *The American Journal of Surgical Pathology*. 2006; 30(2):223-229.
8. Barry TS, Jaffe ES, Sorbara L, et al. Peripheral T-cell lymphomas expressing CD30 and CD15. *The American Journal of Surgical Pathology*. 2003; 27(12):1513-1522.

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Bibliographie

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code : M9-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Mao X, Zhang X, Xue X, et al. Brain tumor stem-like cells identified by neural stem cell marker CD15. Translational Oncology. 2009; 2(4): 247-257.
5. Wong A, Lopategui J, Clancy S, et al. Anaplastic large cell lymphoma associated with a breast implant capsule: a case report and review of the literature. The American Journal of Surgical Pathology. 2008; 32(8):1265-1268.
6. Pellegrini W, Bresciani G, De Zorzi A, et al. MMA monoclonal antibody is a superior anti-CD15 reagent for the diagnosis of classical Hodgkin's lymphoma? Haematologica. 2007; 92(5):708-709.
7. Vassallo J, Lamant L, Brugieres L, et al. ALK positive anaplastic large cell lymphoma mimicking nodular sclerosis Hodgkin's lymphoma: Report of 10 cases. The American Journal of Surgical Pathology. 2006; 30(2):223-229.
8. Barry TS, Jaffe ES, Sorbara L, et al. Peripheral T-cell lymphomas expressing CD30 and CD15. The American Journal of Surgical Pathology. 2003; 27(12):1513-1522.

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Bibliografia

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Mao X, Zhang X, Xue X, et al. Brain tumor stem-like cells identified by neural stem cell marker CD15. Translational Oncology. 2009; 2(4): 247-257.
5. Wong A, Lopategui J, Clancy S, et al. Anaplastic large cell lymphoma associated with a breast implant capsule: a case report and review of the literature. The American Journal of Surgical Pathology. 2008; 32(8):1265-1268.
6. Pellegrini W, Bresciani G, De Zorzi A, et al. MMA monoclonal antibody is a superior anti-CD15 reagent for the diagnosis of classical Hodgkin's lymphoma? Haematologica. 2007; 92(5):708-709.
7. Vassallo J, Lamant L, Brugieres L, et al. ALK positive anaplastic large cell lymphoma mimicking nodular sclerosis Hodgkin's lymphoma: Report of 10 cases. The American Journal of Surgical Pathology. 2006; 30(2):223-229.
8. Barry TS, Jaffe ES, Sorbara L, et al. Peripheral T-cell lymphomas expressing CD30 and CD15. The American Journal of Surgical Pathology. 2003; 27(12):1513-1522.

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Bibliografie

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 28. February 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD und Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Mao X, Zhang X, Xue X, et al. Brain tumor stem-like cells identified by neural stem cell marker CD15. Translational Oncology. 2009; 2(4): 247-257.
5. Wong A, Lopategui J, Clancy S, et al. Anaplastic large cell lymphoma associated with a breast implant capsule: a case report and review of the literature. The American Journal of Surgical Pathology. 2008; 32(8):1265-1268.
6. Pellegrini W, Bresciani G, De Zorzi A, et al. MMA monoclonal antibody is a superior anti-CD15 reagent for the diagnosis of classical Hodgkin's lymphoma? Haematologica. 2007; 92(5):708-709.
7. Vassallo J, Lamant L, Brugieres L, et al. ALK positive anaplastic large cell lymphoma mimicking nodular sclerosis Hodgkin's lymphoma: Report of 10 cases. The American Journal of Surgical Pathology. 2006; 30(2):223-229.
8. Barry TS, Jaffe ES, Sorbara L, et al. Peripheral T-cell lymphomas expressing CD30 and CD15. The American Journal of Surgical Pathology. 2003; 27(12):1513-1522.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code : M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Mao X, Zhang X, Xue X, et al. Brain tumor stem-like cells identified by neural stem cell marker CD15. Translational Oncology. 2009; 2(4): 247-257.
5. Wong A, Lopategui J, Clancy S, et al. Anaplastic large cell lymphoma associated with a breast implant capsule: a case report and review of the literature. The American Journal of Surgical Pathology. 2008; 32(8):1265-1268.
6. Pellegrini W, Bresciani G, De Zorzi A, et al. MMA monoclonal antibody is a superior anti-CD15 reagent for the diagnosis of classical Hodgkin's lymphoma? Haematologica. 2007; 92(5):708-709.
7. Vassallo J, Lamant L, Brugieres L, et al. ALK positive anaplastic large cell lymphoma mimicking nodular sclerosis Hodgkin's lymphoma: Report of 10 cases. The American Journal of Surgical Pathology. 2006; 30(2):223-229.
8. Barry TS, Jaffe ES, Sorbara L, et al. Peripheral T-cell lymphomas expressing CD30 and CD15. The American Journal of Surgical Pathology. 2003; 27(12):1513-1522.

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Kaynakça

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Mao X, Zhang X, Xue X, et al. Brain tumor stem-like cells identified by neural stem cell marker CD15. Translational Oncology. 2009; 2(4): 247-257.
5. Wong A, Lopategui J, Clancy S, et al. Anaplastic large cell lymphoma associated with a breast implant capsule: a case report and review of the literature. The American Journal of Surgical Pathology. 2008; 32(8):1265-1268.
6. Pellegrini W, Bresciani G, De Zorzi A, et al. MMA monoclonal antibody is a superior anti-CD15 reagent for the diagnosis of classical Hodgkin's lymphoma? Haematologica. 2007; 92(5):708-709.
7. Vassallo J, Lamant L, Brugieres L, et al. ALK positive anaplastic large cell lymphoma mimicking nodular sclerositis Hodgkin's lymphoma: Report of 10 cases. The American Journal of Surgical Pathology. 2006; 30(2):223-229.
8. Barry TS, Jaffe ES, Sorbara L, et al. Peripheral T-cell lymphomas expressing CD30 and CD15. The American Journal of Surgical Pathology. 2003; 27(12):1513-1522.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.

4. Mao X, Zhang X, Xue X, et al. Brain tumor stem-like cells identified by neural stem cell marker CD15. *Translational Oncology*. 2009; 2(4): 247-257.
5. Wong A, Lopategui J, Clancy S, et al. Anaplastic large cell lymphoma associated with a breast implant capsule: a case report and review of the literature. *The American Journal of Surgical Pathology*. 2008; 32(8):1265-1268.
6. Pellegrini W, Bresciani G, De Zorzi A, et al. MMA monoclonal antibody is a superior anti-CD15 reagent for the diagnosis of classical Hodgkin's lymphoma? *Haematologica*. 2007; 92(5):708-709.
7. Vassallo J, Lamant L, Brugieres L, et al. ALK positive anaplastic large cell lymphoma mimicking nodular sclerosis Hodgkin's lymphoma: Report of 10 cases. *The American Journal of Surgical Pathology*. 2006; 30(2):223-229.
8. Barry TS, Jaffe ES, Sorbara L, et al. Peripheral T-cell lymphomas expressing CD30 and CD15. *The American Journal of Surgical Pathology*. 2003; 27(12):1513-1522.

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Leica Biosystems Newcastle Ltd
Leica Biosystems Inc
Pty Ltd