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# Novocastra™ Liquid Mouse Monoclonal Antibody

## Beta-Sarcoglycan

### Product Code: NCL-L-b-SARC

#### Intended Use

*For in vitro diagnostic use.*

NCL-L-b-SARC is intended for the qualitative identification by light microscopy of Beta-Sarcoglycan by immunohistochemistry. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

#### Principle of Procedure

Immunohistochemical (IHC) staining techniques allow for the visualization of antigens via the sequential application of a specific antibody to the antigen (primary antibody), a secondary antibody to the primary antibody and an enzyme complex with a chromogenic substrate with interposed washing steps. The enzymatic activation of the chromogen results in a visible reaction product at the antigen site. The specimen may then be counterstained and coverslipped. Results are interpreted using a light microscope and aid in the differential diagnosis of pathophysiological processes, which may or may not be associated with a particular antigen.

#### Clone

βSarc1/5B1

#### Immunogen

Fusion protein RBSG-NT of the human beta-sarcoglycan sequence.

#### Specificity

Human beta-sarcoglycan (43 kD).

#### Reagent Composition

NCL-L-b-SARC is a liquid tissue culture supernatant containing sodium azide as a preservative.

#### Ig Class

IgG1

#### Total Protein Concentration

Total Protein

Refer to vial label for lot specific total protein concentration.

#### Antibody Concentration

Greater than or equal to 26 mg/L as determined by ELISA. Refer to vial label for batch specific Ig concentration.

#### Recommendations On Use

Immunohistochemistry (see **Methodology**) on frozen sections. Suggested dilution: 1:50–1:200 for 60 minutes at 25 °C. This is provided as a guide and users should determine their own optimal working dilutions.

#### Storage and Stability

Store at 2–8 °C. Do not freeze. Return to 2–8 °C immediately after use. Do not use after expiration date indicated on the vial label. Storage conditions other than those specified above must be verified by the user.

#### Specimen Preparation

Freeze specimen tissue blocks in isopentane chilled in liquid nitrogen (see Warnings and Precautions). The specimens do not require further fixation but should be embedded in OCT™ compound (Sakura, Product No. Tissue-Tek 4583).

#### Warnings and Precautions

This reagent has been prepared from the supernatant of cell culture. As it is a biological product, reasonable care should be taken when handling it.

This reagent contains sodium azide. A Material Safety Data Sheet is available upon request or available from [www.LeicaBiosystems.com](http://www.LeicaBiosystems.com) Consult federal, state or local regulations for disposal of any potentially toxic components.

Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions.<sup>1</sup> Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents and specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.

Minimize microbial contamination of reagents or an increase in non-specific staining may occur.

Incubation times or temperatures, other than those specified, may give erroneous results. Any such changes must be validated by the user.

Liquid nitrogen due to its excessively cold temperature causes burns and protective clothing, including gloves and visor, should be used when handling. Use in a well ventilated area.

Isopentane is highly flammable and harmful by ingestion and inhalation. It is also irritating to skin and eyes, and is a narcotic in high concentration.

## Quality Control

Differences in tissue processing and technical procedures in the user's laboratory may produce significant variability in results, necessitating regular performance of in-house controls in addition to the following procedures.

Controls should be fresh autopsy/biopsy/surgical specimens frozen as soon as possible in the same manner as the patient sample(s).

## Positive Tissue Control

Used to indicate correctly prepared tissues and proper staining techniques.

One positive tissue control should be included for each set of test conditions in each staining run.

A tissue with weak positive staining is more suitable than a tissue with strong positive staining for optimal quality control and to detect minor levels of reagent degradation.<sup>2</sup>

Recommended positive control tissue is skeletal muscle.

If the positive tissue control fails to demonstrate positive staining, results with the test specimens should be considered invalid.

## Negative Tissue Control

Should be examined after the positive tissue control to verify the specificity of the labeling of the target antigen by the primary antibody.

Recommended negative control tissue has not been evaluated.

Alternatively, the variety of different cell types present in most tissue sections frequently offers negative control sites, but this should be verified by the user.

Non-specific staining, if present, usually has a diffuse appearance. Sporadic staining of connective tissue may also be observed in sections from excessively formalin-fixed tissues. Use intact cells for interpretation of staining results. Necrotic or degenerated cells often stain non-specifically.<sup>3</sup> False-positive results may be seen due to non-immunological binding of proteins or substrate reaction products. They may also be caused by endogenous enzymes such as pseudoperoxidase (erythrocytes), endogenous peroxidase (cytochrome C), or endogenous biotin (eg. liver, breast, brain, kidney) depending on the type of immunostain used. To differentiate endogenous enzyme activity or non-specific binding of enzymes from specific immunoreactivity, additional patient tissues may be stained exclusively with substrate chromogen or enzyme complexes (avidin-biotin, streptavidin, labeled polymer) and substrate-chromogen, respectively. If specific staining occurs in the negative tissue control, results with the patient specimens should be considered invalid.

## Negative Reagent Control

Use a non-specific negative reagent control in place of the primary antibody with a section of each patient specimen to evaluate non-specific staining and allow better interpretation of specific staining at the antigen site.

## Patient Tissue

Examine patient specimens stained with NCL-L-b-SARC last. Positive staining intensity should be assessed within the context of any non-specific background staining of the negative reagent control. As with any immunohistochemical test, a negative result means that the antigen was not detected, not that the antigen was absent in the cells/tissue assayed. If necessary, use a panel of antibodies to identify false-negative reactions.

## Results Expected

### Normal Tissues

Clone  $\beta$ Sarc1/5B1 detects the beta-sarcoglycan protein in the sarcolemma of skeletal, cardiac and smooth muscle fibers. It also crossreacts with beta-sarcoglycan in sections of muscle from mouse, rat, rabbit, dog, chicken, hamster and pig.

### Abnormal Tissues

Clone  $\beta$ Sarc1/5B1 has been used in immunohistochemical and immunoblotting studies of more than 850 patients to identify a deficiency of the 43 kD dystrophin-associated glycoprotein, beta-sarcoglycan.

**NCL-L-b-SARC is recommended for the identification of human Beta-Sarcoglycan by immunohistochemistry.**

## General Limitations

Immunohistochemistry is a multistep diagnostic process that consists of specialized training in the selection of the appropriate reagents; tissue selection, fixation, and processing; preparation of the IHC slide; and interpretation of the staining results.

Tissue staining is dependent on the handling and processing of the tissue prior to staining. Improper fixation, freezing, thawing, washing, drying, heating, sectioning or contamination with other tissues or fluids may produce artifacts, antibody trapping, or false negative results. Inconsistent results may be due to variations in fixation and embedding methods, or to inherent irregularities within the tissue.<sup>4</sup> Excessive or incomplete counterstaining may compromise proper interpretation of results.

The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Antibodies from Leica Biosystems Newcastle Ltd are for use, as indicated, on either frozen or paraffin-embedded sections with specific fixation requirements. Unexpected antigen expression may occur, especially in neoplasms. The clinical interpretation of any stained tissue section must include morphological analysis and the evaluation of appropriate controls.

## **Bibliography - General**

1. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991; 7(9). Order code M29-P.
2. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. *Progress in Surgical Pathology*. 6:1–15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
3. Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. *Laboratory Medicine*. 1983; 14:767.
4. Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. *American Journal of Clinical Pathology*. 1980; 73:626.
5. Anderson LVB. Multiplex Western blot analysis of the muscular dystrophy proteins. Chapter 22, p369–386, in *Muscular Dystrophy: Methods and Protocols* (number 43 in the *Methods in Molecular Medicine* series), edited by Bushby KMD & Anderson LVB. 2001. Humana Press: Totowa, New Jersey.
6. Pogue R, Anderson LVB, Pyle A, et al. Strategy for mutation analysis in the autosomal recessive limb-girdle muscular dystrophies. *Neuromuscular Disorders*. 2001; 11(1):80–87.
7. Auranen M, Rapola J, Pihko H, et al. Muscle membrane-skeleton protein changes and histopathological characterization of muscle-eye-brain disease. *Neuromuscular Disorders*. 2000; 10(1):16–23.
8. Bornemann A and Anderson LVB. Diagnostic protein expression in human muscle biopsies. *Brain Pathology*. 2000; 10:193–214.
9. Vainzof M, Moreira ES, Canovas M, et al. Partial alpha-sarcoglycan deficiency with retention of the dystrophin-glycoprotein complex in a LGMD2D family. *Muscle Nerve*. 2000; 23(6):984–988.
10. Anderson LVB and Davison K. Multiplex Western blotting system for the analysis of muscular dystrophy proteins. *American Journal of Pathology*. 1999; 154(4):1017–1022.
11. O'Brien DP, Johnson GC, Liu LA, et al. Laminin alpha2 (merosin) –deficient muscular dystrophy and demyelinating neuropathy in two cats. *Journal of the Neurological Sciences* 2001; 189:37–43.
12. Sheriffs IN, Rampling D and Smith VV. Paraffin wax embedded muscle is suitable for the diagnosis of muscular dystrophy. *Journal of Clinical Pathology*. 2001; 54(7):517–520.
13. Ueda H, Ueda K, Baba T, et al. Delta and gamma-sarcoglycan localization in the sarcoplasmic reticulum of skeletal muscle. *Journal of Histochemistry and Cytochemistry*. 2001; 49:529–538.

## **Amendments to Previous Issue**

Not Applicable.

## **Date of Issue**

17 October 2018

# Immunohistochemistry methodology for using Novocastra™ antibodies on frozen muscle tissue.

## Reagents required but not supplied

1. Standard solvents used in immunohistochemistry.
2. 50 mM Tris-buffered saline (TBS) pH 7.6.
3. Antibody diluent - normal serum optimally diluted in TBS.
4. Normal serum from the species in which the secondary antibody is raised.
5. Secondary peroxidase-conjugated antibody - use as recommended by manufacturer.
6. 3,3' Diaminobenzidine tetrahydrochloride (DAB) - prepare and use as recommended by manufacturer.
7. Mounting medium - use as recommended by manufacturer.

## Equipment required but not supplied

1. Incubator set to 25 °C.
2. General immunohistochemistry laboratory equipment.
3. Electric fan for air drying slides.

## Antigen retrieval solutions (see Recommendations on Use)

Not applicable to frozen sections.

## Methodology

Prior to undertaking this methodology, users must be trained in immunohistochemical techniques.

Users should determine optimal dilutions for antibodies. Unless indicated, all steps are performed at 25 °C.

1. Cut and mount 4–10 µm sections on slides coated with a suitable tissue adhesive and air dry for at least one hour.
2. Incubate sections with optimally diluted primary antibody (see Recommendations on Use).
3. Wash in TBS buffer for 2 x 5 minutes with gentle rocking.
4. Incubate sections in appropriate peroxidase-conjugated secondary antibody.
5. Wash in TBS buffer for 2 x 5 minutes with gentle rocking.
6. Incubate the slides in DAB.
7. Rinse the slides in clean water.
8. Dehydrate, clear and mount sections.

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### **Literatur - Allgemein**

1. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991;7(9). Order code M29-P.
2. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. *Progress in Surgical Pathology*. 6:1–15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
3. Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. *Laboratory Medicine*. 1983; 14:767.
4. Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. *American Journal of Clinical Pathology*. 1980; 73:626.
5. Anderson LVB. Multiplex Western blot analysis of the muscular dystrophy proteins. Chapter 22, p369–386, in *Muscular Dystrophy: Methods and Protocols* (number 43 in the *Methods in Molecular Medicine* series), edited by Bushby KMD & Anderson LVB. 2001. Humana Press: Totowa, New Jersey.
6. Pogue R, Anderson LVB, Pyle A, et al. Strategy for mutation analysis in the autosomal recessive limb-girdle muscular dystrophies. *Neuromuscular Disorders*. 2001; 11(1):80–87.
7. Auranen M, Rapola J, Pihko H, et al. Muscle membrane-skeleton protein changes and histopathological characterization of muscle-eye-brain disease. *Neuromuscular Disorders*. 2000; 10(1):16–23.
8. Bornemann A and Anderson LVB. Diagnostic protein expression in human muscle biopsies. *Brain Pathology*. 2000; 10:193–214.
9. Vainzof M, Moreira ES, Canovas M, et al. Partial alpha-sarcoglycan deficiency with retention of the dystrophin-glycoprotein complex in a LGMD2D family. *Muscle Nerve*. 2000; 23(6):984–988.
10. Anderson LVB and Davison K. Multiplex Western blotting system for the analysis of muscular dystrophy proteins. *American Journal of Pathology*. 1999; 154(4):1017–1022.
11. O'Brien DP, Johnson GC, Liu LA, et al. Laminin alpha2 (merosin) –deficient muscular dystrophy and demyelinating neuropathy in two cats. *Journal of the Neurological Sciences* 2001; 189:37–43.
12. Sheriffs IN, Rampling D and Smith VV. Paraffin wax embedded muscle is suitable for the diagnosis of muscular dystrophy. *Journal of Clinical Pathology*. 2001; 54(7):517–520.
13. Ueda H, Ueda K, Baba T, et al. Delta and gamma-sarcoglycan localization in the sarcoplasmic reticulum of skeletal muscle. *Journal of Histochemistry and Cytochemistry*. 2001; 49:529–538.

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### **Bibliografia - Geral**

1. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991; 7(9). Order code M29-P.
2. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. *Progress in Surgical Pathology*. 6:1–15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F, Field & Wood, Inc., Philadelphia.
3. Nadjji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. *Laboratory Medicine*. 1983; 14:767.
4. Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. *American Journal of Clinical Pathology*. 1980; 73:626.
5. Anderson LVB. Multiplex Western blot analysis of the muscular dystrophy proteins. Chapter 22, p369–386, in *Muscular Dystrophy: Methods and Protocols* (number 43 in the *Methods in Molecular Medicine* series), edited by Bushby KMD & Anderson LVB. 2001. Humana Press: Totowa, New Jersey.
6. Pogue R, Anderson LVB, Pyle A, et al. Strategy for mutation analysis in the autosomal recessive limb-girdle muscular dystrophies. *Neuromuscular Disorders*. 2001; 11(1):80–87.
7. Auranen M, Rapola J, Pihko H, et al. Muscle membrane-skeleton protein changes and histopathological characterization of muscle-eye-brain disease. *Neuromuscular Disorders*. 2000; 10(1):16–23.
8. Bornemann A and Anderson LVB. Diagnostic protein expression in human muscle biopsies. *Brain Pathology*. 2000; 10:193–214.
9. Vainzof M, Moreira ES, Canovas M, et al. Partial alpha-sarcoglycan deficiency with retention of the dystrophin-glycoprotein complex in a LGMD2D family. *Muscle Nerve*. 2000; 23(6):984–988.
10. Anderson LVB and Davison K. Multiplex Western blotting system for the analysis of muscular dystrophy proteins. *American Journal of Pathology*. 1999; 154(4):1017–1022.
11. O'Brien DP, Johnson GC, Liu LA, et al. Laminin alpha2 (merosin) –deficient muscular dystrophy and demyelinating neuropathy in two cats. *Journal of the Neurological Sciences* 2001; 189:37–43.
12. Sheriffs IN, Rampling D and Smith VV. Paraffin wax embedded muscle is suitable for the diagnosis of muscular dystrophy. *Journal of Clinical Pathology*. 2001; 54(7):517–520.
13. Ueda H, Ueda K, Baba T, et al. Delta and gamma-sarcoglycan localization in the sarcoplasmic reticulum of skeletal muscle. *Journal of Histochemistry and Cytochemistry*. 2001; 49:529–538.

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## Bibliografi - Allmän

1. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991; 7(9). Order code M29-P.
2. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. *Progress in Surgical Pathology*. 6:1–15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
3. Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and its pitfalls. *Laboratory Medicine*. 1983; 14:767.
4. Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. *American Journal of Clinical Pathology*. 1980; 73:626.
5. Anderson LVB. Multiplex Western blot analysis of the muscular dystrophy proteins. Chapter 22, p369–386, in *Muscular Dystrophy: Methods and Protocols* (number 43 in the *Methods in Molecular Medicine* series), edited by Bushby KMD & Anderson LVB. 2001. Humana Press: Totowa, New Jersey.
6. Pogue R, Anderson LVB, Pyle A, et al. Strategy for mutation analysis in the autosomal recessive limb-girdle muscular dystrophies. *Neuromuscular Disorders*. 2001; 11(1):80–87.
7. Auranen M, Rapola J, Pihko H, et al. Muscle membrane-skeleton protein changes and histopathological characterization of muscle-eye-brain disease. *Neuromuscular Disorders*. 2000; 10(1):16–23.
8. Bornemann A and Anderson LVB. Diagnostic protein expression in human muscle biopsies. *Brain Pathology*. 2000; 10:193–214.
9. Vainzof M, Moreira ES, Canovas M, et al. Partial alpha-sarcoglycan deficiency with retention of the dystrophin-glycoprotein complex in a LGMD2D family. *Muscle Nerve*. 2000; 23(6):984–988.
10. Anderson LVB and Davison K. Multiplex Western blotting system for the analysis of muscular dystrophy proteins. *American Journal of Pathology*. 1999; 154(4):1017–1022.
11. O'Brien DP, Johnson GC, Liu LA, et al. Laminin alpha2 (merosin) –deficient muscular dystrophy and demyelinating neuropathy in two cats. *Journal of the Neurological Sciences* 2001; 189:37–43.
12. Sheriffs IN, Rampling D and Smith VV. Paraffin wax embedded muscle is suitable for the diagnosis of muscular dystrophy. *Journal of Clinical Pathology*. 2001; 54(7):517–520.
13. Ueda H, Ueda K, Baba T, et al. Delta and gamma-sarcoglycan localization in the sarcoplasmic reticulum of skeletal muscle. *Journal of Histochemistry and Cytochemistry*. 2001; 49:529–538.

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1. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991; 7(9). Order code M29-P.
2. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. Progress in Surgical Pathology. 6:1–15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
3. Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. Laboratory Medicine. 1983; 14:767.
4. Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. American Journal of Clinical Pathology. 1980; 73:626.
5. Anderson LVB. Multiplex Western blot analysis of the muscular dystrophy proteins. Chapter 22, p369–386, in Muscular Dystrophy: Methods and Protocols (number 43 in the Methods in Molecular Medicine series), edited by Bushby KMD & Anderson LVB. 2001. Humana Press: Totowa, New Jersey.
6. Pogue R, Anderson LVB, Pyle A, et al. Strategy for mutation analysis in the autosomal recessive limb-girdle muscular dystrophies. Neuromuscular Disorders. 2001; 11(1):80–87.
7. Auranen M, Rapola J, Pihko H, et al. Muscle membrane-skeleton protein changes and histopathological characterization of muscle-eye-brain disease. Neuromuscular Disorders. 2000; 10(1):16–23.
8. Bornemann A and Anderson LVB. Diagnostic protein expression in human muscle biopsies. Brain Pathology. 2000; 10:193–214.
9. Vainzof M, Moreira ES, Canovas M, et al. Partial alpha-sarcoglycan deficiency with retention of the dystrophin-glycoprotein complex in a LGMD2D family. Muscle Nerve. 2000; 23(6):984–988.
10. Anderson LVB and Davison K. Multiplex Western blotting system for the analysis of muscular dystrophy proteins. American Journal of Pathology. 1999; 154(4):1017–1022.
11. O'Brien DP, Johnson GC, Liu LA, et al. Laminin alpha2 (merosin) –deficient muscular dystrophy and demyelinating neuropathy in two cats. Journal of the Neurological Sciences 2001; 189:37–43.
12. Sheriffs IN, Rampling D and Smith VV. Paraffin wax embedded muscle is suitable for the diagnosis of muscular dystrophy. Journal of Clinical Pathology. 2001; 54(7):517–520.
13. Ueda H, Ueda K, Baba T, et al. Delta and gamma-sarcoglycan localization in the sarcoplasmic reticulum of skeletal muscle. Journal of Histochemistry and Cytochemistry. 2001; 49:529–538.

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Leica Biosystems Inc
Pty Ltd