

NEW COLORS AVAILABLE

## PERMASLIDE and PERMASLIDE Plus Thermal Slides

Designed and validated specifically for the new HistoCore PERMA S Slide Printer, the PERMASLIDE portfolio addition provides consistent, high quality, easy to read imprints for use with thermal print technology.

- New colors available for easy color-coding system to help reduce errors and standardize workflow.
- Adhesive clipped or non-clipped corner
- Non-adhesive clipped or non-clipped corner
- Optimized for thermal printing



Advancing Cancer Diagnostics  
Improving Lives

**Leica**  
BIO SYSTEMS

Product Number	Description
3800449	White Non-Clipped PERMASLIDE
3800455	White Non-Clipped PERMASLIDE Plus Adhesive
3800449CL	White Clipped PERMASLIDE
3800451CL	Pink Clipped PERMASLIDE
3800453CL	Blue Clipped PERMASLIDE
3800454CL	Yellow Clipped PERMASLIDE
3800447CL	Orange Clipped PERMASLIDE
3800448CL	Tan Clipped PERMASLIDE
3800455CL	White Clipped PERMASLIDE Plus Adhesive
3800456CL	Pink Clipped PERMASLIDE Plus Adhesive
3800458CL	Blue Clipped PERMASLIDE Plus Adhesive
3800459CL	Yellow Clipped PERMASLIDE Plus Adhesive
3800444CL	Orange Clipped PERMASLIDE Plus Adhesive
3800443CL	Tan Clipped PERMASLIDE Plus Adhesive

Each Product Quantity: 1 case (1,000 slides)



Contact your local Leica Biosystems representative to learn more about our PERMASLIDE and PERMASLIDE Plus Thermal Slides