

N o w A v a i l a b l e

Mega Cassettes, Double Depth

Expand Your Tissue Processing Capabilities
with the Recently Released Mega Cassettes,
Double Depth

FREEDOM TO PROCESS SPECIMENS OF ALL SIZES

- These cassettes offer double the depth of conventional cassettes.
- Allows for easier embedding and processing of whole organs, ossified and larger tissues (such as eyes, teeth and bones) that do not fit in standard cassettes.
- Perfect for research and specialty clinical applications.
- The cassette dimensions are 40 x 26 x 13 mm and are available in 6 colors in quantities of 1,000.
- Can be used with the Mega cassettes base mold made of stainless steel. These reusable base molds allow easy block removal after cooling. Available in boxes of 10.



| Product Number | Description |
|----------------|--------------------------------------|
| 38VSP59040 | Mega Cassettes, Double Depth, White |
| 38VSP59041 | Mega Cassettes, Double Depth, Orange |
| 38VSP59042 | Mega Cassettes, Double Depth, Blue |
| 38VSP59044 | Mega Cassettes, Double Depth, Yellow |
| 38VSP59045 | Mega Cassettes, Double Depth, Green |
| 38VSP59046 | Mega Cassettes, Double Depth, Red |
| 3803085E | 32 x 23 x 14 mm Metal Base Mold |

Contact your local Leica Biosystems representative to learn more about our newly available Mega Cassettes.

LeicaBiosystems.com

Copyright © 2019 Leica Biosystems Richmond, Inc. All Rights Reserved. LEICA and the Leica logo are registered trademarks of Leica Microsystems IR GmbH.

181630 Rev B 09/2019

Leica

B I O S Y S T E M S

Advancing Cancer Diagnostics
Improving Lives