

# Potassium Chloride (0.075M)

---

**Cat. No:** KBI-90052 (100 ml)

**Store at:** 2 - 8°C

## Product Description

Potassium Chloride 0.075 M is a hypotonic solution, utilized for the preparation of blood lymphocyte chromosomes.

## Instructions for Use

Work according to own established protocols.

## Precaution and Disclaimer

1. For *in vitro* diagnostic use. It is not intended for therapeutic use.
2. In case of emergencies check SDS for medical advice.
3. Use of Potassium Chloride does not guarantee the successful outcome of any diagnostic testing.
4. Do not use Potassium Chloride beyond the expiration date indicated on the product label.

## Storage and Stability

Potassium Chloride should be kept at 2-8°C, protected from light



**Kreatech  
Biotechnology B.V.**  
Vlierweg 20  
1032 LG Amsterdam  
The Netherlands