

Sodium Citrate Solution (0.8%)

Cat. No: KBI-90054 (500 ml)

Store at: 2 - 8°C

Product Description

Sodium Citrate Solution (0.8%) is a hypotonic solution, utilized for the preparation of blood lymphocyte chromosomes.

Precaution and Disclaimer

1. For *in vitro* diagnostic use. It is not intended for therapeutic use.
2. In case of emergencies check SDS for medical advice.
3. Use of Sodium Citrate Solution does not guarantee the successful outcome of any diagnostic testing.
4. Do not use Sodium Citrate Solution beyond the expiration date indicated on the product label.

Storage and Stability

Sodium Citrate Solution should be kept at 2-8°C, protected from light

Instructions for Use

Work according to own established protocols.



**Kreatech
Biotechnology B.V.**
Vlierweg 20
1032 LG Amsterdam
The Netherlands