

# Novocastra™ Liquid Mouse Monoclonal Antibody CD10

**Product Code: NCL-L-CD10-270**

<b>Intended Use</b>	FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.
<b>Specificity</b>	Human CD10 molecule, also known as common acute lymphocytic leukemia antigen (CALLA).
<b>Clone</b>	56C6
<b>Ig Class</b>	IgG1
<b>Antigen Used for Immunizations</b>	Prokaryotic recombinant fusion protein corresponding to the external domain of the CD10 glycoprotein.
<b>Hybridoma Partner</b>	Mouse myeloma (p3-NS1-Ag4-1).
<b>Preparation</b>	Liquid tissue culture supernatant containing sodium azide. Volume as indicated on vial label.
<b>Effective on Frozen Tissue</b>	Yes. Acetone fixation recommended.
<b>Effective on Paraffin Wax Embedded Tissue</b>	Yes
<b>Recommendations on Use</b>	Immunohistochemistry on paraffin sections. <b>Heat Induced Epitope Retrieval (HIER):</b> Please follow the instructions for use in Novocastra Epitope Retrieval Solution pH 6. <b>Suggested dilution:</b> 1:100 for 30 minutes at 25 °C. This is provided as a guide and users should determine their own optimal working dilutions. <b>Visualization:</b> Please follow the instructions for use in the Novolink™ Polymer Detection Systems. For further product information or support, contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems web site, <a href="http://www.LeicaBiosystems.com">www.LeicaBiosystems.com</a> <u>The performance of this antibody should be validated when utilized with other manual staining systems or automated platforms.</u> <b>Western Blotting:</b> Typical working dilution 1:50–1:100.
<b>Positive Controls</b>	Immunohistochemistry: Tonsil. Western Blotting: Placenta.
<b>Staining Pattern</b>	Membrane
<b>Storage and Stability</b>	Store at 2–8 °C. Do not freeze. Return to 2–8 °C immediately after use. Do not use after expiration date indicated on the vial label. Storage conditions other than those specified above must be verified by the user.
<b>Warnings and Precautions</b>	This reagent has been prepared from the supernatant of cell culture. As it is a biological product, reasonable care should be taken when handling it. This reagent contains sodium azide. A Material Safety Data Sheet is available upon request or available from <a href="http://www.LeicaBiosystems.com">www.LeicaBiosystems.com</a>



**General Overview**

CD10 is a 100 kD cell surface metalloendopeptidase called neprilysin which inactivates a variety of biologically active peptides. It was initially identified as the common acute lymphoblastic leukemia antigen (CALLA) and considered to be tumor-specific. Subsequent studies, however, have shown that CD10 is expressed on the surface of a wide variety of normal and neoplastic cells. CD10 has been identified on the surface of normal early lymphoid progenitor cells, immature B cells within adult bone marrow and germinal centre B cells within lymphoid tissue. It is also expressed in various non-lymphoid cells and tissues, such as breast myoepithelial cells, bile canaliculi, fibroblasts, with especially high expression on the brush border of kidney and gut epithelial cells.

**General References**

Ordi J, Romagosa C, Tavassoli FA et al. *American Journal of Surgical Pathology* 2003 27(2), 178–186.  
Kristiansen G, Schluns K, Yongwei Y et al. *Annals of Cellular Pathology* 2002 24(1), 41–46.  
Sumathi V P and McCluggage WG. *Journal of Clinical Pathology* 2002 55, 391–392.  
Dunphy CH, Polsky JM, Lance Evans H et al. *Leuk Lymphoma* 2001 41(5–6), 585–592.  
Eshoa C, Perkins S, Kampalath et al. *American Journal of Clinical Pathology* 2001 115(6), 862–867.  
Tajima Y, Shimoda T, Nakanishi Y et al. *Oncology* 2001 61(3), 212–220.  
Watson P, Wood KM, Lodge A et al. *Histopathology* 2000 36(2), 145–150.  
Bavikatty NR, Ross CW, Finn WG et al. *Human Pathology* 2000 31(9), 1051–1054.  
Chen CC, Raikow RB, Sonmez-Alpan E et al. *Appl Immunohistochem Mol Morphol* 2000 8(1), 1–11.  
Chu P and Arber DA. *American Journal of Clinical Pathology* 2000 113(3), 374–382.  
Chu PG, Chang KL, Weiss LM et al. *Appl Immunohistochem Mol Morphol* 2000 8(4), 257–262.  
Conde-Sterling DA, Aguilera NS, Nandedkar MA et al. *Arch Pathol Lab Med* 2000 124(5), 704–708.  
Endoh Y, Tamura G, Motoyama T et al. *Human Pathology* 1999 30(7), 826–832.  
Kaufmann O, Flath B, Späth-Schwalbe E et al. *American Journal of Clinical Pathology* 1999 111(1), 117–122.  
McIntosh GG, Lodge AJ, Watson P et al. *American Journal of Pathology* 1999 154(1), 77–82.  
Millar E K, Waldron S, Spencer A et al. *Journal of Clinical Pathology* 1999 52, 849–850.