

Novocastra™ Predilute Rabbit Polyclonal Antibody *Helicobacter pylori*

Product Code: IHC5994-A

Analyte Specific Reagent

Intended Use	This reagent is an analyte specific reagent. Analytical and performance characteristics are not established.		
Reagents Provided	Helicobacter pylori is a predilute rabbit polyclonal antibody presented in Tris buffer pH 7.7 with 10mg/ml carrier protein and 0.1% ProClin® 300. Total volume = 7 mL.		
Specificity	Reactive with heat stable, somatic antigens of the whole <i>Helicobacter pylori</i> organism.		
Warnings and Precautions	<p>Helicobacter pylori Contains a mixture of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC No. 247-500-7] and 2-methyl-2H-isothiazol-3-one [EC No. 220-239-6] (3:1) (<1%). Signal words: Warning.</p>	<p>H317: May cause an allergic skin reaction.</p>	<p>P261: Avoid breathing dust/fumes/gas/mist/vapours/spray. P272: Contaminated work clothing should not be allowed out of the workplace. P280: Wear protective gloves/protective clothing/eye protection/face protection. P302+352: IF ON SKIN: Wash with plenty of water/soap. P321: Specific treatment (see instructions on this label). P333+313: If skin irritation or rash occurs: Get medical attention. P362+364: Take off contaminated clothing and wash it before reuse. P501: Dispose of contents/container to hazardous or special waste collection point.</p>
Statement of Quality	Each lot of reagent has been quality controlled by immunohistochemistry.		
Storage and Stability	Store at 2–8 °C. This product is stable under these conditions up to the expiry date indicated on the vial label.		
References	National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A, 1991; 7(9). Order code M29-P.		
Date of Issue	23 February 2016.		

