

Novocastra Predilute Rabbit Polyclonal Antibody *Helicobacter pylori*

Product Code: IHC5994-A

Analyte Specific Reagent

Intended Use This reagent is an analyte specific reagent. Analytical and performance characteristics are not established.

Reagents Provided *Helicobacter pylori* is a predilute rabbit polyclonal antibody presented in Tris buffer pH 7.7 with 10mg/ml carrier protein and 0.1% ProClin® 300.
Total volume = 7 mL.

Specificity Reactive with heat stable, somatic antigens of the whole *Helicobacter pylori* organism.

Warnings and Precautions

Helicobacter pylori
Contains 0.035%
2-methylisothiazol-3(2H)-
one.
GHS07: Exclamation
mark.
Signal words: Warning.

H317: May cause an
allergic skin reaction.

P261: Avoid breathing dust/fumes/gas/mist/
vapours/spray.
P272: Contaminated work clothing should not be
allowed out of the workplace.
P280: Wear protective gloves/protective clothing/
eye protection/face protection.
P302+P352: IF ON SKIN: Wash with plenty of
water/soap.
P321: Specific treatment (see instructions on
this label).
P333+P313: If skin irritation or rash occurs: Get
medical attention.
P362+P364: Take off contaminated clothing and
wash it before reuse.
P501: Dispose of contents/container to hazardous
or special waste collection point.

Statement of Quality

Each lot of reagent has been quality controlled by immunohistochemistry.

Storage and Stability

Store at 2–8 °C. This product is stable under these conditions up to the expiry date indicated on the vial label.

References

National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A, 1991; 7(9). Order code M29-P.

Date of Issue

25 February 2022

