

Novocastra™ Liquid Mouse Monoclonal Antibody Immunoglobulin A

Product Code: NCL-L-IgA

Intended Use	FOR RESEARCH USE ONLY.
Specificity	Alpha chain of human IgA
Clone	N1CLA
Ig Class	IgG1
Antigen Used for Immunizations	Prokaryotic recombinant protein corresponding to 129 amino acids of the heavy chain constant region of the human immunoglobulin A molecule.
Hybridoma Partner	Mouse myeloma (p3-NS1-Ag4.1)
Preparation	Liquid tissue culture supernatant containing 15 mM sodium azide. Volume as indicated on vial label.
Effective on Frozen Tissue	Not evaluated.
Effective on Paraffin Wax Embedded Tissue	Yes (using heat induced epitope retrieval with citrate-based buffer, pH 6.0: see overleaf)
Recommendations on Use	Immunohistochemistry: Typical working dilution 1:300. Heat induced epitope retrieval technique using Citrate-based buffer, pH 6.0. 30 minutes primary antibody incubation at 25 °C. Polymer detection recommended.
Positive Controls	Immunohistochemistry: Tonsil
Staining Pattern	Cytoplasmic
Storage and Stability	Store liquid antibody at 4 °C. Under these conditions, there is no significant loss in product performance up to the expiry date indicated on the vial label. Prepare working dilutions on the day of use.
General Overview	IgA is a member of the antibody class of the immunoglobulin superfamily. There are several classes and subclasses (isotypes) of antibody, the antibody isotype being defined by the immunoglobulin heavy chain present in the molecule. The basic structure of an immunoglobulin molecule consists of two identical heavy chains ($\gamma, \mu, \alpha, \zeta, \epsilon$) and two identical light chains, either kappa or lambda. IgA contains the α chain and may be present in a serum or secretory form. In serum, 90% of IgA is monomeric, while in its secretory form it is the main immunoglobulin found in secretions including tears, saliva, intestinal and bronchial mucus, sweat, colostrum, and secretions from the prostate and respiratory epithelia, where it has the job of defending exposed external surfaces of the body against attack from microorganisms. Secretory IgA is synthesized locally by plasma cells and dimerized intracellularly with a cysteine rich J-chain.
General References	Fagarasan S and Honjo T. Current Opinion in Immunology. 16: 277–283 (2004). Pilette C, Ouadrhiri Y, Godding V et al. European Respiratory Journal. 18: 571–588 (2001).



Instructions for Use

Heat Induced Epitope Retrieval Combined With Polymer Detection For Immunohistochemical Demonstration On Paraffin Sections

1. Cut and mount sections on slides coated with a suitable tissue adhesive.
2. Deparaffinize sections and rehydrate to distilled water.
3. Place sections in 0.5% hydrogen peroxide/methanol for 10 minutes (or use other appropriate endogenous peroxidase blocking procedure). Wash sections in tap water.
4. Heat 1500 mL of the recommended epitope retrieval solution (Citrate based pH 6.0 - Epitope Retrieval Solution unless otherwise indicated overleaf) in a stainless steel pressure cooker until boiling. Cover but do not lock lid.
5. Position slides into metal staining racks (do not place slides close together as uneven staining may occur) and lower into pressure cooker ensuring slides are completely immersed in epitope retrieval solution. Lock lid.
6. When the pressure cooker reaches operating temperature and pressure (after about 5 minutes) start a timer for 1 minute (unless otherwise indicated on the data sheet).
7. When the timer rings, remove pressure cooker from heat source and run under cold water with lid on. **DO NOT OPEN LID UNTIL THE INDICATORS SHOW THAT PRESSURE HAS BEEN RELEASED.** Open lid, remove slides and place immediately into a bath of tap water.
8. Wash sections once using fresh Tris-Buffered Saline (TBS, pH 7.6) buffer for 5 minutes.
9. Place sections in diluted normal serum (eg NCL-G-SERUM) for 10 minutes.
10. Incubate sections with primary antibody.
11. Wash twice, each time using fresh TBS buffer for 5 minutes.
12. For visualization of the bound primary antibody, follow instructions supplied with the Polymer Detection System.
13. Counterstain with hematoxylin (if required), dehydrate and mount.

** (In most applications, Phosphate Buffered Saline, pH 7.6, can be used instead of TBS, pH 7.6).*

Safety Note

To ensure the correct and safe use of your pressure cooker, PLEASE READ THE MANUFACTURER'S INSTRUCTIONS.