



Surgipath® SelectTech®

Hematoxylin and Eosin Staining System

Surgipath® has undertaken the task of standardizing the reagents and stains used in the hematoxylin and eosin staining procedure. The reagents and stains that make up this new system must pass vigorous analytical chemistry standards and performance tests. The built-in Quality Assurance parameters assure consistent results. The reagents and stains in this system were designed to complement each other. The result is a crisp, well delineated hematoxylin and eosin stain.

- Excellent cost per slide ratio.
- Two year shelf life.

Surgipath® SelecTech®

Hematoxylin and Eosin Staining System

- Surgipath® SelecTech Hematoxylin 560™ a progressive hematoxylin that stains darker the longer it is immersed (up to 5 minutes).
Can be used for routine, IHC counterstaining or for frozen section staining.
- Surgipath SelecTech Hematoxylin 560 MX™ a progressive hematoxylin that produces clearly delineated nuclear chromatin with darker results. Great for those who prefer darker nuclei.
- Surgipath Define® provides sharp definition of cytoplasmic details and increases cellular clarity. Mucin and ground substance background staining are removed with Define.
- Surgipath Define MX-aq™ an aqueous differentiating reagent designed to work with the Surgipath SelecTech Hematoxylin 560 MX and Surgipath SelecTech Hematoxylin 560, but very effective with other hematoxylin.
- Surgipath SelecTech Eosin Phloxine 515™ produces a more reddish coloration in the cytoplasm, especially in connective tissue components.
- Surgipath SelecTech Eosin 515 LT™ produces a lighter staining cytoplasm with a more pastel appearance.
- Surgipath Blue Buffer 8™ at a pH of 8, assures proper bluing of nuclei.

Product Name	Packaging	Product Code
Surgipath® SelecTech Hematoxylin 560™	500 mL, each	3801570
Surgipath® SelecTech Hematoxylin 560™	500 mL, 4/case	3801571
Surgipath® SelecTech Hematoxylin 560 MX™	500 mL, 4/case	3801575
Surgipath® SelecTech Eosin Phloxine 515™	500 mL, 4/case	3801606
Surgipath® SelecTech Eosin 515 LT™	500 mL, 4/case	3801619
Surgipath® SelecTech Eosin Trichrome™	500 mL, each	3801610
Surgipath® SelecTech Eosin Trichrome™	500 mL, 4/case	3801611
Surgipath® SelecTech Alcoholic Eosin Y 515™	500 mL, each	3801615
Surgipath® SelecTech Alcoholic Eosin Y 515™	500 mL, 4/case	3801616
Surgipath® Define®	500 mL, each	3803590
Surgipath® Define®	500 mL, 4/case	3803591
Surgipath® Define MX-aq™	500 mL, each	3803595
Surgipath® Define MX-aq™	500 mL, 4/case	3803596
Surgipath® Blue Buffer 8™	500 mL, each	3802915
Surgipath® Blue Buffer 8™	500 mL, 4/case	3802916

For More Information,

Call Your Local Sales Representative

Leica Biosystems Richmond, Inc.
 5205 Route 12, Richmond, Illinois 60071
 Tel. 815-678-2000
 CS 800-225-3035 Fax 815-678-2216

Leica Biosystems Canada

83 Terracon Place, Winnipeg, Manitoba R2J 4B3
 Tel. 204-632-7150 800-665-7425 Fax 204-632-7093