



BIOSYSTEMS

ACD RNA ISH Probe (concentrate) HPV-11

Product Code: RS7860

Analyte Specific Reagent

Intended Use	This reagent is an analyte specific reagent. Analytical and performance characteristics are not established.
Reagent Provided	HPV-11 Probe is an HPV 11 oligonucleotide probe at >2000 ng/mL supplied in TE buffer. Minimum volume = 80 µL.
Specificity	Human papillomavirus type 11 E6, E7 mRNA Transcript.
Species	<i>Human papillomavirus type 11</i>
Target Gene	E6, E7
Alternate Names	N/A
Entrez Gene ID	N/A
Accession Number	FR872717.1
Target Region	149-761
Warnings and Precautions	For additional information on toxicity please refer to the Safety Data Sheet (SDS), which is available upon request; alternatively visit www.LeicaBiosystems.com
Statement of Quality	Each lot of reagent has been quality controlled by probe hybridization to a complementary plasmid sequence.
Storage and Stability	Store at 2–8 °C. The product is stable under these conditions up to the expiry date indicated on the vial label.
Date of Issue	15 October 2020

